



Features:

- + Titanium tack material, good biocompatibility and MRI friendly.
- + Low profile titanium mesh interface (no cap).
- + Matte finish shaft, no reflective light in procedure.
- + Two working lengths available to meet the demands of laparoscopic and open surgery.
- + Ergonomic pistol handle design, for comfort and easy manipulation.

Code No.	Item	Spec.(mm)	Units/ Box
MFT-O-01	Hernia Stapler for open surgery	5 x 110	1
MFT-L-02	Hernia Stapler for laparoscopic surgery	5 x 350	1