


Hydrophilic Guide Wire



Features

- Innovative triple anti-drop design at Hydrophilic tip, with no risk of drop-off.
- Anti-twist innerNiti core wire, offering an excellent twisting and pushing force.
- Super smooth PTFE zebra coating, easier to pass through the working channel.
- Yellow/Black coating, easier to track the guide wire and obvious under X-ray.
- Unique design of tapering tip, highly increase the success rate of surgery.

Product code	O.D.(mm)	Length(mm)	Characteristics	Units/Box
CIT-GW-80-150-C1	0.80 (0.032")	1500		5
CIT-GW-80-260-C1	0.80 (0.032")	2600		5
CIT-GW-80-450-C1	0.80 (0.032")	4500	C1: Hydrophilic Guide	5
CIT-GW-88-150-C1	0.88 (0.035")	1500	Wire.	5
CIT-GW-88-260-C1	0.88 (0.035")	2600	Bicolor spiral; Super	5
CIT-GW-88-450-C1	0.88 (0.035")	4500	smooth soft head; Straight	5
CIT-GW-63-150-C1	0.63(0.025")	1500	shape and angled shape	5
CIT-GW-63-260-C1	0.63(0.025")	2600		5
CIT-GW-63-450-C1	0.63(0.025")	4500		5