






# Implantable Ligating Clips



## Features:

- + Three sizes to ligate vessel structures from 3mm to 16mm.
- + Integrated teeth to ensure greater stability and prevent slippage on the vessel.
- + Distinct tactile feedback confirms jaw seating and secure vessel placement.
- + Clips delivered in color-coded cartridge for easy identification.
- + No interference with CT, MRI or X-ray diagnostics.

Code No.	Specification	Color of Cartridge	Cartridges/ Box
LC-01	XL	Yellow 	14
LC-02	L	Purple 	14
LC-03	ML	Green 	14

## Implantable Ligating Clips

Application	Surgery	Surgical Method and Clip Size Selection
General Surgery	Colectomy: Left, Right or Middle Colon Resection	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
General Surgery	Anterior Rectal Resection	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
General Surgery	Appendicectomy	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
General Surgery	Cholecystectomy	Endoscopic Surgery: ML, L Open Surgery: ML, L
General Surgery	Gastrectomy	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
General Surgery	Fundoplication	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
General Surgery	Splenectomy	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
Urology	Prostatectomy	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
Urology	Adrenalectomy	Endoscopic Surgery: ML, L Open Surgery: ML, L
Urology	Nephrectomy	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
Gynecology	Hysterectomy	Endoscopic Surgery: ML, L, XL Open Surgery: ML, L, XL
Thoracic Surgery	Pneumonectomy, Pulmonary lobectomy, Segmentectomy	Endoscopic Surgery: XL Open Surgery: ML, L, XL
Vascular Surgery	Venous Tumor Resection	Open Surgery: ML, L Endoscopic Surgery: XL
Cardiovascular Surgery	Coronary Artery Bypass Graft	Open Surgery: ML, L, XL
Cardiovascular Surgery	Abdominal Arterial Aneurysm Resection	Endoscopic Surgery: XL Open Surgery: ML, L, XL
Otolaryngology	Thyroidectomy	Open Surgery: ML, L