


Injection Needle



Features

- Innovative design of metallic tip, efficiently avoid damage to catheter and working channel.
- Rigid and super-lubricant catheter, easier to pass through the working channel.
- Unique needle section design, with suitable sharpness and easier for puncture. Ergonomic handle design, easier to operate.

Product code	O .D.(mm) (D)	Working Length (mm) (L)	Needle size (diameter/ length) (d/l)	Endoscope Channel (mm)	Units/Box
CIT-IN-DLdl-A	1.8	1200,1800	19G, 21G,23G, 25G, 4/5/6mm	≥2.0	10
CIT-IN-DLdl-A	2.3	1800,2300	19G, 4/5/6mm,21G, 4/5/6mm 23G, 4/5/6mm,25G, 4/5/6mm	≥2.8	10