

Irrigation Tube



Green Cap
For Asia-Pacific Baxter Sterile Water Bottle Use



Transparent cap
For B.Braun and Fresenius Sterile Water Bottle Use

Features

- With back-flow valve, provide irrigation via sterile water during GI endoscopic procedures.
- Always connect with brand new water jet connector, the irrigation tube can be used for 24 hours. Discard daily.
- Universal bottle cap compatible for most water bottles on the market, like B.Braun, Fresenius and Baxter.
- Two kinds of pump tubing, compatible with most irrigation units: Olympus® AFU-100, Olympus® OFP, Olympus® OFP-2, Medivators EGP-100, ERBE EIP2 or Endotechnik., etc.

Product Code	Description	Units/Box
CIT-IT-005	Irrigation Tube (For Olympus AFU-100, Olympus OFP, Medivators EGP-100, ERBE EIP2 or Endotechnik Irrigation Units), diameter of pump tube:6.4mm. Compatible with B.Braun and Fresenius sterile water bottles.	20
CIT-IT-005C	Irrigation Tube (For Olympus OFP and OFP-2 Irrigation Units), diameter of pump tube:9.7mm. Compatible with B.Braun and Fresenius sterile water bottles.	20
CIT-IT-010	Irrigation Tube (For Olympus AFU-100, Olympus OFP, Medivators EGP-100, ERBE EIP2 or Endotechnik Irrigation Units), diameter of pump tube:6.4mm. Compatible with Asia-Pacific Baxter sterile water bottles.	20
CIT-IT-010C	Irrigation Tube (For Olympus OFP and OFP-2 Irrigation Units), diameter of pump tube:9.7mm. Compatible with Asia-Pacific Baxter sterile water bottles.	20