



Intended use:

Prevent the spread of germs in hospitals and in geriatric nursing home; for medical personnel, patients and visitors use in the ward.

Features:

- + Comply to ANSI/AAMI PB70:2012
- + Coated PP+ PE material provide great water proof performance and and air permeability. A highly effective barrier against fluids, germs and light chemical splashes.
- + Ultrasonic welding technology used on the sleeves of isolation gowns.
- + Classification levels of barrier performance: Level 1 to 3 available.
- + Size of isolation gowns: from S to 4XL.
- + Color: Blue/White/Yellow.

Code No.	Name	Size	Color	Raw Material
CIT-GLF-L1	Isolation Gown – Level 1	S/M/L/XL/2XL/3XL/4XL	Blue/White/Yellow	35gsm PP+PE
CIT-GLF-L2	Isolation Gown – Level 2	S/M/L/XL/2XL/3XL/4XL	Blue/White/Yellow	40gsm PP+PE
CIT-GLF-L2	Isolation Gown – Level 3	S/M/L/XL/2XL/3XL/4XL	Blue/White/Yellow	40gsm PP+PE
CIT-GLF-28	Isolation Gown – Unrated	S/M/L/XL/2XL/3XL/4XL	Blue/White/Yellow	28gsm PP
CIT-GLF-35	Isolation Gown – Unrated	S/M/L/XL/2XL/3XL/4XL	Blue/White/Yellow	35gsm PP

Packaging:

Pcs/Bag	Bags/CTN	Pcs/CTN	Carton Size(cm)
25pcs	4bags	100pcs	52x32x42

• Products can be customized. Please contact us for any further information.

