


J Shape Fistula Occluding Tracheal Stent



Application: J Shape Fistula-Occluding Tracheal stent is used for tracheal and bronchial stenosis caused by various malignant lesions.

Performance Characteristics:

- J-design fitting anatomic structure
- Self-expanding with high radial force
- Excellent positional stability
- Resistant and elastic covering
- Fully covered stents available
- Release with high radiopacity
- Golden radiopaque markers for clearer visibility



Specifications:

Product Code	Tracheal/ Bronchial Dia (mm)	Tracheal/ Bronchial Length(mm)	Product Illustration	Introducer
NTJA2-FT1-16030	16/12	30/30		18Fr*600mm
NTJA2-FT2-18040	18/14	40/30		
NTJA2-FT3-20050	20/14	50/30		
NTJA2-FT4-16030	16/12	30/30		24Fr*700mm
NTJA2-FT5-18040	18/14	40/30		
NTJA2-FT6-20050	20/14	50/30		