

Hemorrhoids Stapler & STARR Stapler with Transparent Housing



Effective Volume of the

Housing Increases 35%!



CE0197

Hemorrhoids Stapler & STARR Stapler with Transparent Housing



Ref. No.	Outside Diameter (mm)	Cutting Diameter (mm)	No. of Staples (Piece)	Staple Height (mm)
PAH32-TF	32	22.2	30	3.8
PAH34-TF	34	24.2	32	3.8

STARR Stapler

Ref. No.	Outside Diameter (mm)	Cutting Diameter (mm)	No. of Staples (Piece)	Staple Height (mm)
PAH32-TF-ST-4.5	32	22.2	30	4.5
PAH34-TF-ST-4.5	34	24.2	32	4.5

Features:

- + The large volume staple housing allows the resection of more tissue when clinically required.
- + Transparent housing provides better visibility during the procedure.
- + Stainless steel made shaft ensures the stability during the clamping and firing process.
- + Clear cutting voice reminds the surgeon the firing is completely done.
- + Detachable or fixed anvils for surgeon's preference.
- + Butterfly shape adjusting knob, easier for manipulation.
- + Safety Mechanism helps to avoid the misfiring.
- + Four Suture hole design for better access and easier operation.
- + Transparent anoscope with scale for better observation.