


Polymer Clip Appliers






Right Applier Jaw

Keep jaw angle 90° to the cartridge

Features:




- + Detachable and non-detachable, two options offered.
- + Both styles with flushable shaft for ease of cleaning and maintenance.
- + Matte finish shaft, avoid the reflect light in procedure.
- + Color coded size for easy identification.
- + Compatible with corresponding sizes of ligating clips.
- + Different angle jaws available, facilitate firm clip application and better access in the most difficult procedures and position.
- + Wide product range used in both Laparoscopic and open surgery.

Code No.	Specification	Clip size	Rotator Color	Units/ Box
LC-LA-01	10mm x 330mm	XL		1
LC-LA-02	10mm x 330mm	L		1
LC-LA-03	5mm x 330mm	ML		1






20° Polymer Clip Appliers



Code No.	Specification	Clip size	Rotator Color	Units/ Box
LC-LA-01-20	10mm x 330mm, 20°	XL		1
LC-LA-02-20	10mm x 330mm, 20°	L		1
LC-LA-03-20	5mm x 330mm, 20°	ML		1

45° Polymer Clip Appliers



Code No.	Specification	Clip size	Rotator Color	Units/ Box
LC-LA-01-45	10mm x 330mm, 45°	XL		1
LC-LA-02-45	10mm x 330mm, 45°	L		1
LC-LA-03-45	5mm x 330mm, 45°	ML		1

Detachable in 2 parts Polymer Clip Appliers



Code No.	Specification	Clip size	Rotator Color	Units/ Box
LC-LA-01-2	10mm x 330mm, 20°	XL		1
LC-LA-02-2	10mm x 330mm, 20°	L		1
LC-LA-03-2	5mm x 330mm, 20°	ML		1