


Polymer Clip Appliers



Features:

- + Detachable and non-detachable, two options offered.
- + Both styles with flushable shaft for ease of cleaning and maintenance.
- + Matte finish shaft, avoid the reflect light in procedure.
- + Color coded size for easy identification.
- + Compatible with corresponding sizes of ligating clips.
- + Different angle jaws available, facilitate firm clip application and better access in the most difficult procedures and position.
- + Wide product range used in both Laparoscopic and open surgery.

Product Name	Specification	Clip size	Rotator Color	Ref. Code No.
Polymer Clip Applier	10mm x 330mm	XL		LC-LA-01
	10mm x 330mm	L		LC-LA-02
	5mm x 330mm	ML		LC-LA-03
Polymer Clip Applier Detachable in 2 parts	10mm x 330mm	XL		LC-LA-01-2
	10mm x 330mm	L		LC-LA-02-2
	5mm x 330mm	ML		LC-LA-03-2

