




Reusable Bulldog Forceps and Bulldog Clips



Straight Type



Curved Type

Intended Use:

Reusable Bulldog Forceps and Reusable Bulldog Clips are designed for the temporary occlusion of vessel during resection procedure.



Specification of Reusable Bulldog Forceps

Ref.Code	Diameter (d)	Working Length (L)	Angle (α)
BD-10380	10mm	380mm	≥15°
BD-12380	12.5mm	380mm	≥15°

Specification of Reusable Bulldog Clips

Ref.Code	Type	Clip Length	With Teeth	Closing Force
BDC-25ZT3(5)	Straight	25mm	Yes	3N/5N
BDC-45ZT3(5)	Straight	45mm	Yes	3N/5N
BDC-65ZT3(5)	Straight	65mm	Yes	3N/5N
BDC-25WT3(5)	Curved	25mm	Yes	3N/5N
BDC-45WT3(5)	Curved	45mm	Yes	3N/5N
BDC-65WT3(5)	Curved	65mm	Yes	3N/5N
BDC-25Z3(5)	Straight	25mm	No	3N/5N
BDC-45Z3(5)	Straight	45mm	No	3N/5N
BDC-65Z3(5)	Straight	65mm	No	3N/5N
BDC-25W3(5)	Curved	25mm	No	3N/5N
BDC-45W3(5)	Curved	45mm	No	3N/5N
BDC-65W3(5)	Curved	65mm	No	3N/5N