


Reusable Nathanson Liver Retractor Set



Specification of the Liver Retractor

Size	Length(mm)	Width(mm)
XL	235	95
L	195	75
M	165	65
S	135	55

Intended Use:

Reusable Nathanson Liver Retractor Set is intended to be used to retract the left lobe of the liver during surgery in the upper gastrointestinal area.

Ref.Code	Composition of the Set
N-LIV-SET	Holder*1
	Mechanical Arm Assembly: Post*1, Cross Bar*1, Hand wheel handle*1, Big Screw*1, Big Clamp*1, Small Screw*1, Plum blossom handle*1
	XL Size Liver Retractor*1
	L Size Liver Retractor*1
	M Size Liver Retractor*1
S Size Liver Retractor*1	