


# Segmented Esophageal Stent



**Application:** Malignant esophageal tumors in esophageal middle segment, distal segment and some cases on the upper segment

## Performance Characteristics:

- Segmented modular stent body allowing independent movements between each module to maximize conformance to curved to tortuous anatomy
- Segmented design minimizes migration
- Comparably mild radial force and low axial force
- Elastic covering to resist tumor inner growth
- Extraction threads for easy removal and repositioning
- Delivery system with markers for better position accuracy



## Specifications:

Product Code	Dia. A (mm)	Total Length B (mm)	Product Illustration	Introducer
NTJB3-AB- FC	18/20	70/80/85/100/120/135/150		24Fr*700mm OTW
NTJB3-AB- ARFC				

Remark: FC means full covered; ARFC means full covered and anti reflux.