

Anti-Skid Shoe cover



Intended use:

Wear shoe covers to provide a barrier against possible exposure to airborne organisms or contact with a contaminated environment. Use shoe covers for patients with hemorrhagic disease. Shoe covers should be worn as part of Full Barrier Precautions.

Features:

- + Made of non-woven fabric, dirt and dust and other particles cannot stick to the surface and be ground into the weave.
- + Anti-skid design protect you against spills and biohazards, while keeping you against slippage on a wet floor.
- + Good elasticity and rebound effect, durable.

Code No.	Size(cm)	Color	Raw Material
CIT-FHXT-15	40×15	Blue	PP Non-woven fabric, 30±5gsm, Anti-skid
CIT-FHXT-17	40×17	Blue	PP Non-woven fabric, 30±5gsm, Anti-skid

Packaging:

Pcs/Bag	Bags/CTN	Pcs/CTN	Carton Size(cm)
150pcs	14bags	2100pcs	72x28x42

• Products can be customized. Please contact us for any further information.

