



Intended use:

Surgical gowns are worn by medical staff during operations, experiments, and research. They are used to isolate infections caused by germs and other factors, avoid secondary infections to patients, and prevent liquids such as alcohol and blood from infecting medical staff. Mainly used for medical operation and patient treatment; epidemic prevention and inspection in public places; disinfection in virus-contaminated areas.

Features:

- + 35gsm or 45gsm SMS non-woven fabric raw material, provide good air permeability, softness and adsorbability.
- + Comply to EN13795-1:2019 and ANSI/AAMI PB70-2012.
- + Knitted cuff design provide comfortable wearing experience for long time operation process.
- + Neck with Velcro design, ensure no lose of the gown.
- + 4 belts design. Tie all the belts to provide completely covering to the medical staff, avoid potential infection.
- + Both sterile and non-sterile surgical gowns are available.

Code No.	Name	Size	Color	Raw Material
CIT-SSY-S	Surgical Gown	S/M/L/XL/2XL/3XL/4XL	Blue	45gsm SMS
CIT-SSY-R	Reinforced Surgical Gown	S/M/L/XL/2XL/3XL/4XL	Blue	35gsm SMS

Packaging:

	Pcs/Pouch	Pouches/CTN	Pcs/CTN	Carton Size(cm)
Sterile (with 2 hand towels)	1pc	25Pouches	25pcs	47x40x36

	Pcs/Bag	Bags/CTN	Pcs/CTN	Carton Size(cm)
Non sterile	10pcs	6bags	60pcs	60x25x48

• Products can be customized. Please contact us for any further information.

