


Titanium Clip Appliers



Features:

- + Detachable and non-detachable, two options offered.
- + Both styles with flushable shaft for ease of cleaning and maintenance.
- + Matte finish shaft, avoid the reflect light in procedure.
- + Compatible with corresponding sizes of titanium clips.
- + Single action, double action and side load 90° tips available.

Product Name	Type	Specification	Clip Size	Code No.	Units/ Box
Titanium Clip Applier for Laparoscopic Surgery	single action	10mm *330mm	L	CIT-RCA-01	1
	single action	10mm *330mm	ML	CIT-RCA-02	1
	single action	5mm *330mm	M	CIT-RCA-03	1
	single action	5mm *330mm	S	CIT-RCA-04	1
	side load 90°	10mm *330mm	ML	CIT-RCA-05	1
	double action	10mm *330mm	ML	CIT-RCA-06	1
	single action and detachable in 2 parts	10mm *330mm	L	CIT-RCA-07	1
	single action and detachable in 2 parts	10mm *330mm	ML	CIT-RCA-08	1
	side load 90° and detachable in 2 parts	10mm *330mm	ML	CIT-RCA-09	1
	double action and detachable in 2 parts	10mm *330mm	L	CIT-RCA-10	1
	double action and detachable in 2 parts	10mm *330mm	ML	CIT-RCA-11	1



Product Name	Specification	Clip Size	Code No.	Units/ Box
Titanium Clip Applier for Open Surgery	267mm	L	CIT-RCA-12	1
	267mm	ML	CIT-RCA-13	1
	191mm	M	CIT-RCA-14	1
	146mm	S	CIT-RCA-15	1

