


# Tracheal Stent

**Application:** Tracheal stent is used for tracheal stenosis caused by various malignant lesions.




Vertical



### Performance Characteristics:

- Self expanding with high radial force
- Excellent positional stability
- Resistant and elastic covering
- Fully covered stents available
- Release with high radiopacity
- Golden radiopaque markers for clearer visibility

### Stent Shape Y



- 1. Vertical cylinder structure
- 17. Tracheal two ball ends

### Specifications:

Product Code	Dia.A (mm)	Total Length B (mm)	Stent Shape Y	Introducer
NTJA2-ABY- FC	10/12/14/16/18/20/22	20/30/40/50/60/70/80	1/17	18Fr*600mm OTW

Remark: FC means full covered.