


# UES Esophageal Stent



**Application: Upper Esophageal Sphincter Cancer**

### Performance Characteristics:

- Proximal special design to alleviate foreign body sensation
- Elastic covering to resist tumor inner growth
- Golden radiopaque markers for clearer visibility
- Distal and proximal deployment
- Delivery system with markers for better position accuracy
- Customized stent is recommended



### Specifications:

Product Code	Dia. A (mm)	Total Length B (mm)	Product Illustration	Introducer
NTJB1-AB- HC	18/20	60/80/100/120/140/150		24Fr*700mm OTW
	22/24			18Fr*650mm Proximal Release
NTJB1-AB- FC	18/20	60/80/100/120/140/150		24Fr*700mm OTW
	22/24			18Fr*650mm Proximal Release

Remark: HC means half covered; FC means full covered.