

Water Jet Connector



CIT-WJC-241



CIT-WJC-251



CIT-WJC-261



CIT-WJC-271

Features

- Single use, eliminate the work burden of reprocessing.
- No lubricant on the valve to interfere the procedure or obstruct the view.
- Water jet connector(CIT-WJC-241) with metal threading design, to minimize potential to cross thread and avoid leakage.
- Valve comes packaged in sterile pouch for ease use.

Product Code	Description	Units/Box
CIT-WJC-241	Water jet connector for Olympus Endoscope(Metal Insert)	100
CIT-WJC-251	Water jet connector for Olympus Endoscope(All Plastic)	100
CIT-WJC-261	Water jet connector for Pentax 90 Series Endoscope	100
CIT-WJC-271	Water jet connector for Fujifilm Endoscope	100