

Disposable Endoscope Valves System For Fujifilm Endoscopes 500/600 Series



Air/Water valve



Suction valve





Biopsy valve



Water jet connector

Features

- Single use, eliminate the work burden of reprocessing.
- Valves come packaged in sterile kits for ease use.
- Great blister packing design enables the valves to drop out of box easily.
- Compatible with Fujifilm® Endoscopes 500/600 series.

Product Picture	Ref. Code	Description	Units/Box
	CIT-VK-016	3-piece Valve kit (Air/Water Valve, Suction Valve, Biopsy Valve)	25
	CIT-VK-017	4-piece Valve kit (Air/Water Valve, Suction Valve, Biopsy Valve, Water Jet Connector)	25