

# Disposable Gas Water Valve

## For Olympus Endoscopes

New



### Intended Use

Disposable Gas Water Valve is intended to used in combination with GI endoscope to supply water and gas of carbon dioxide(CO<sub>2</sub>).

### Features

- Disposable Gas Water Valve is specialized designed for CO<sub>2</sub> insufflation systems.
- Reduce the emission of CO<sub>2</sub> into indoor environment during the procedure.
- Lower the carbon footprint of endoscopy and ensure the endoscopy procedures more environmentally friendly.
- Compatible with Olympus® Endoscopes 140/240/260/290/1200/1500 Series.
- Single use, eliminate the work burden of reprocessing.
- No lubricant on the valves to interfere the procedure or obstruct the view.

Product Picture	Ref. Code	Description	Units/Box
	CIT-GWV-A2	Disposable Gas Water Valve	100
	CIT-VK-018	3-piece Valve Kit (Biopsy Valve, Gas Water Valve, Suction Valve)	25
	CIT-VK-019	4-piece Valve Kit (Biopsy Valve, Gas Water Valve, Suction Valve, Water Jet Connector)	25
	CIT-VK-020	5-piece Valve Kit (Biopsy Valve, Gas Water Valve, Suction Valve, Water Jet Connector, AW Channel Cleaning Adapter)	25

