

# Sterile Tray Liner



## Intended Use

Sterile Tray Liner is designed to be used with endoscope transport tray, to provide full coverage of the tray.

It's composed of:

- 1 clear base liner,
- 1 green liner indicates clean endoscope,
- 1 red liner indicates contaminated endoscope.



## Features

- Sterile tray liner greatly helps to reduce the potential risk of cross-infection.
- Thick base liner with strong anti-burst strength, provide safe environment for endoscope.
- BIOHAZARD and CLEAN markings on corresponding red liner and green liner, provides clear distinguish between contaminated and clean endoscope.
- Transparent color of red and green liners, helps to identify the inside packed endoscope easily.
- Individually packed in blister pouch, reduce 50% plastic waste compared to PE pouch.

Product Code	Description	Units/Box
EZE-LINER	Sterile, compatible with CLEANASCOPE style tray	200